-Scratches or dents can result in breakage.
-For instruments produced by another manufacturer, reference the manufacturer's instructions for use.
-Care should be taken to remove debris, tissue or bone fragments that may collect on the instrument.
LIMITATIONS ON REPROCESSING:
-End of useful life is generally determined by wear or damage in surgical use.
-Carefully inspect instruments between uses to verify proper functioning and configuration.
-Damaged instruments must be replaced to prevent potential patient injury such as loss of metal fragments into the surgical site.

POSTERIOR SPINE TRUSS SYSTEM (PSTS) INTERBODY FUSION DEVICE

INSTRUCTIONS FOR USE (ENGLISH)

PLEASE READ CAREFULLY

CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALES BY OR ON THE ORDER OF A PHYSICIAN.

-For patients with, or suspected with, Creutzfeldt-Jakob Disease (CJD), variant CJD or other transmissible spongiform encephalopathy (TSE) and related infections, it is recommended to treat the patient using single-use instruments.

INTENDED USE:
The Posterior Spine Truss System (PSTS) Interbody Fusion Device is designed to provide mechanical support to the lumbar spine while biologic tissue takes place.

DEVICE DESCRIPTION:
The device is an open architecture truss design mathematically formulated to provide structural support with open space throughout the implant for bone growth and fusion. The 4WEB additive manufacturing process provides a hierarchical surface roughness. The implant is made from Ti6Al4V alloy.

The device is available in a variety of sizes and lordotic angles to accommodate the patient's anatomy.

INDICATIONS FOR USE:

The Posterior Spine Truss System (PSTS) Interbody Fusion Device is indicated for use in skeletal maturity patients with Degenerative Disc Disease (DDD) at one or two contiguous levels from L2-S1. DDD is defined as discogenic back pain with degeneration of the disc confirmed by patient history and radiographic studies. Patients should have received 6 months of non-operative treatment prior to treatment with the devices. The device must be implanted in a supplemental fusion and must be used with autograft and/or allogeneic bone graft composed of cancellous and/or corticocancellous bone graft. These DDD patients may also have up to Grade 1 spondylolisthesis or retrolisthesis at the involved level(s).

CONTRAINDICATIONS:

The PSTS Interbody Fusion Device should not be implanted in patients with:
-An active infection at the operative site or other active systemic infections
-Tumor involvement at the operative site
-Prior fusion at the level(s) to be treated
-Known sensitivity to the material

WARNINGS AND PRECAUTIONS:

-The following warnings have not been established in the following conditions:
-Obesity
-Smoking
-Three or more levels to be fused
-Symptomatic cardiac disease
-Pregnancy
-Previous fusion attempts at the involved level(s)
-Spondylolisthesis or retrolisthesis greater than Grade 1
-Significant loss of bone stock as well as osteoporosis or osteomalacia
-Conditions requiring chronic corticosteroid use
-Active drug abuse

-Place the PSTS Interbody Fusion Device for single use only.
-Bending or fracture of the implants or instruments can occur if not handled properly.
-All sizes and instrumentation are provided non-sterile and must be steam sterilized prior to use. The 4WEB Spine Truss System Implants are intended to be used together with the 4WEB Spine Truss System instrumentation during spinal fusion interventions.

MAINTENANCE INSPECTION AND FUNCTIONAL TESTING:
-Visually inspect all instruments to ensure no damage and wear.
-Ensure there are no cracked handles and shafts are secure in handles.
-Ensure long instruments are free of any bending and distortion.
-Ensure instrument tips are free of defects or burrs.
-Ensure complex instruments with moving parts function appropriately.

CLEANING AND DISINFECTION:
-Inspect all instruments to verify that all visible debris is removed during cleaning and prior to sterilization. If debris is still visible after cleaning, repeat the cleaning process.
-Use the instrument only for the intended purpose.

WRAPPING TECHNIQUE:
-Use instrument trays to contain instruments that are provided in sets.
-Biologic or Chemical Indicators (Bi or Chi) used for monitoring the performance of sterilization processes should be placed in the middle racks within wrapped trays. They should be tested according to the Bi or Chi manufacturer's directions.
-Double wrap instruments in accordance with local procedures, using standard wrapping techniques such as those described in ANSI/AAMI ST79.

OPERATIVE PRECAUTIONS:
The surgeon is to be thoroughly familiar with the PSTS Interbody Fusion Device, methods of application, indications and surgical technique. Correct positioning of the PSTS Interbody Fusion Device relative to the vertebrae should be checked intraoperatively with x-ray. The size (and more particularly the height) of the PSTS Interbody Fusion Device must be chosen on the basis of the patient's anatomy and desired correction.

Each PSTS Interbody Fusion Device is to be filled with autologous and/or allograft bone to promote bone fusion (See Surgical Technique Manual for complete details). The implants are for single-implant use only. An explanted implant must never be re-implanted. Stresses and fracture, even though not noticeable by visual inspection, may have been created during initial implantation.

POTENTIAL ADVERSE EVENTS:
Potential adverse events related to surgery in general, spine surgery specifically or the device. These may include, but are not limited to the following:
-Adverse events related to any surgery: reactions to anesthesia, the anesthetic or other medications; bleeding; infection; wound-healing problems; tissue damage; atelectasis; pneumonia; hematoma; seroma; wound dehiscence or incisional hernia; urologic problems; embolism; anemia; colitis; thrombophlebitis; heart attack; stroke; death.
-Adverse events related specifically to spine surgery, dual tear and CSF leak; nerve damage leading to radiculopathy, myelopathy, paraparesis, paresthesia or paralysis; meningitis; vertebral body damage or fracture; ligament damage; fractured sacrum; or retrograde ejaculation.
-Adverse events related to the device: implant crack or fracture, failure to achieve fusion, implant migration, dislodgement, or metal sensitivity to a foreign body, including possible tumor formation. Additional surgery may be necessary for implant removal, repositioning or replacement. Additional stabilization at the implanted area may be necessary if another disc level may be necessary if non-union or anatomic change at an adjacent level develops.

POSTOPERATIVE CARE:
The physician's postoperative directions and warnings to the patient and the corresponding patient compliance are extremely important.
-The use and limitations of the device must be given to the patient. The patient must be warned that loosening and / or breakage of the device(s) are complications which may occur as result of early or excessive weight-bearing, muscular activity or sudden jolts or shock to the spine.
-The patient must be advised not to smoke or consume alcohol during period of the bone fusion process.
-The patient must be advised of the inability to bend at the point of spinal fusion and taught to compensate for this permanent physical restriction in body motion.
-The patient must be informed and confirmed by roentgenographic examination, if a non-union develops or if the components loosen, migrate, and / or break, the device must be revised and / or removed immediately before serious injury occurs.
-Any revised devices must not be used in any other surgical procedure.

MR SAFETY INFORMATION:
-Non-clinical testing has demonstrated the PSTS Interbody Fusion Device is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:
-Static magnetic field of 3 T or 1.5 T
-Spatial gradient of 100 Gauss/cm or equivalent (19 T/m)
-Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, non-clinical testing results indicate the PSTS Interbody Fusion Device provides a maximum temperature increase of no more than 5.0°C after 15 minutes of continuous scanning.
-In non-clinical testing, the image artifact caused by the device extends approximately 28 mm from the PSTS Interbody Fusion Device when imaged with a gradient echo pulse sequence and a 3 T MR system.

PACKAGING:
PSTS Interbody Fusion Device is provided sterile and is clearly labeled as such in an unopened sterile package provided by 4WEB. The contents are considered sterile unless the package is damaged, opened, or the expiration date on the device label has passed. The integrity of the packaging should be checked to ensure that the sterility of the contents is not compromised.

Implants supplied sterilized from 4WEB must not be re-sterilized.
PSTS ancillary instrumentation is provided non-sterile and is supplied in a tray that is to be steam sterilization prior to use in the operating room. In the case of instrument restock, individual items will be sent in a plastic bag, labeled for that device.

Storage conditions must maintain the integrity of the implants, associated ancillary instruments and their respective packaging. The condition of all implants and instruments must be checked before use. Damaged products must not be used and should be returned to 4WEB.

INFORMATION FOR CLEANING AND STERILIZATION OF SURGICAL INSTRUMENTS:
CAUTION. THESE INSTRUCTIONS DO NOT APPLY TO SINGLE-USE DEVICES.
The instruments used to implant the PSTS Interbody Fusion Device do not have an indefinite functional life. All reusable instruments are subjected to repeated stresses related to bone contact, impact, friction, loading, cleaning and sterilization processes. Instruments should be carefully inspected before use to ensure that they are fully functional.

All ancillary instruments of the PSTS Interbody Fusion Device are delivered non-sterile and therefore, must be decontaminated, cleaned and sterilized prior to surgical use. Decontamination reduces the population of microorganisms and facilitates the subsequent cleaning stage. Strict compliance with the instructions for use pertaining to decontamination and cleaning is mandatory, particularly the concentration and exposure time requirements. Thorough rinsing with water must be conducted following decontamination and cleaning.

-Decontamination: Each hospital must use their own validated decontamination procedures.
-Cleaning: Wash instruments in a LANCER type or equivalent ISO 15883-1 and ISO 15883-2 compliant washing machine with the appropriate cleaning products, rinse, and dry. Any product which may alter the material is prohibited, i.e. bleach, formalin, hypochlorite solutions, saline solution, etc.

The reprocessing instructions have been validated as being capable of preparing reusable 4WEB instruments. It is the responsibility of the reprocessor to ensure that the reprocessing is actually performed according to the manufacturer's instructions and personnel to achieve the desired result. This normally requires validation and routine monitoring of the process. Any deviation by the reprocessor from these instructions should be evaluated for effectiveness and potential adverse consequences.

WARNINGS:
-Follow the instructions and warnings issued by the suppliers of any cleaning and disinfection agents and equipment used.
-Do not exceed 140°C (284°F) during reprocessing steps.
-Highly alkaline conditions can damage products with aluminum parts.
-Acid exposure to hypochlorite solutions, as these will promote corrosion.

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