




POSTERIOR SPINE *TRUSS* SYSTEM™

4WEB's *TRUSS* Implant Technology™ redefined the standard of care as the first FDA-cleared 3D printed interbody fusion device with an Advanced Structural Design. The PSTS continues the legacy of delivering innovation, strength, and performance to every procedure.



THERAPEUTIC BY DESIGN™



- ▶ Truss Implant Technology™ provides a synergistic relationship between the hierarchical surface roughness and the macro architecture.¹
- ▶ After surface features initiate osteogenic activity, the structural mechanics of the truss implant provide a mode of action that converts physiologic load into therapeutic strain at the surgical site, amplifying osteogenic activity in local osteo-biologic material.²
- ▶ This therapeutic effect provides ongoing support while biologic fusion takes place.³

Open Architecture Design allows for bone growth throughout the entire implant construct³

Snowshoe Interface distributes load across the bone contact surface, reducing the risk of subsidence³

Static and articulating inserter options to aid proper implant placement

Bi-convex anatomical design to increase endplate contact



¹ Rowe, et al., SMISS, Annual Forum '19, p.52


² Lee et al. Adv. Healthcare Mater. 2025, 2501856

³ Data on file

PRODUCT CATALOG



VISIT
PRODUCT
PAGE



PSTS IMPLANT SPECIFICATIONS

CATALOG #	FOOTPRINT (W, L, H)	LORDOSIS
STRAIGHT		
PSTS-SM0006-22-SP	9 x 22 x 6mm	0°
PSTS-SM0008-22-SP	9 x 22 x 8mm	0°
PSTS-SM0010-22-SP	9 x 22 x 10mm	0°
PSTS-SM0012-22-SP	9 x 22 x 12mm	0°
PSTS-SM0014-22-SP	9 x 22 x 14mm	0°
PSTS-SM0608-22-SP	9 x 22 x 8mm	6°
PSTS-SM0609-22-SP	9 x 22 x 9mm	6°
PSTS-SM0610-22-SP	9 x 22 x 10mm	6°
PSTS-SM0611-22-SP	9 x 22 x 11mm	6°
PSTS-SM0612-22-SP	9 x 22 x 12mm	6°
PSTS-SM0613-22-SP	9 x 22 x 13mm	6°
PSTS-SM0614-22-SP	9 x 22 x 14mm	6°
PSTS-SM1210-22-SP	9 x 22 x 10mm	12°
PSTS-SM1211-22-SP	9 x 22 x 11mm	12°
PSTS-SM1212-22-SP	9 x 22 x 12mm	12°
PSTS-SM1213-22-SP	9 x 22 x 13mm	12°
PSTS-SM1214-22-SP	9 x 22 x 14mm	12°
PSTS-SM0006-26-SP	9 x 26 x 6mm	0°
PSTS-SM0008-26-SP	9 x 26 x 8mm	0°
PSTS-SM0010-26-SP	9 x 26 x 10mm	0°
PSTS-SM0012-26-SP	9 x 26 x 12mm	0°
PSTS-SM0014-26-SP	9 x 26 x 14mm	0°
PSTS-SM0608-26-SP	9 x 26 x 8mm	6°
PSTS-SM0609-26-SP	9 x 26 x 9mm	6°
PSTS-SM0610-26-SP	9 x 26 x 10mm	6°
PSTS-SM0611-26-SP	9 x 26 x 11mm	6°
PSTS-SM0612-26-SP	9 x 26 x 12mm	6°
PSTS-SM0613-26-SP	9 x 26 x 13mm	6°
PSTS-SM0614-26-SP	9 x 26 x 14mm	6°
PSTS-SM1210-26-SP	9 x 26 x 10mm	12°
PSTS-SM1211-26-SP	9 x 26 x 11mm	12°
PSTS-SM1212-26-SP	9 x 26 x 12mm	12°
PSTS-SM1213-26-SP	9 x 26 x 13mm	12°
PSTS-SM1214-26-SP	9 x 26 x 14mm	12°
PSTS-MD0608-26-SP	11 x 26 x 8mm	6°
PSTS-MD0609-26-SP	11 x 26 x 9mm	6°
PSTS-MD0610-26-SP	11 x 26 x 10mm	6°
PSTS-MD0611-26-SP	11 x 26 x 11mm	6°
PSTS-MD0612-26-SP	11 x 26 x 12mm	6°
PSTS-MD0613-26-SP	11 x 26 x 13mm	6°
PSTS-MD0614-26-SP	11 x 26 x 14mm	6°
PSTS-MD1210-26-SP	11 x 26 x 10mm	12°
PSTS-MD1211-26-SP	11 x 26 x 11mm	12°

CATALOG #	FOOTPRINT (W, L, H)	LORDOSIS
PSTS-MD1212-26-SP	11 x 26 x 12mm	12°
PSTS-MD1213-26-SP	11 x 26 x 13mm	12°
PSTS-MD1214-26-SP	11 x 26 x 14mm	12°
PSTS-MD0608-30-SP	11 x 30 x 8mm	6°
PSTS-MD0609-30-SP	11 x 30 x 9mm	6°
PSTS-MD0610-30-SP	11 x 30 x 10mm	6°
PSTS-MD0611-30-SP	11 x 30 x 11mm	6°
PSTS-MD0612-30-SP	11 x 30 x 12mm	6°
PSTS-MD0613-30-SP	11 x 30 x 13mm	6°
PSTS-MD0614-30-SP	11 x 30 x 14mm	6°
PSTS-MD1210-30-SP	11 x 30 x 10mm	12°
PSTS-MD1211-30-SP	11 x 30 x 11mm	12°
PSTS-MD1212-30-SP	11 x 30 x 12mm	12°
PSTS-MD1213-30-SP	11 x 30 x 13mm	12°
PSTS-MD1214-30-SP	11 x 30 x 14mm	12°
CURVED		
TSTS-MD0007-27-SP	10 x 27 x 7mm	0°
TSTS-MD0009-27-SP	10 x 27 x 9mm	0°
TSTS-MD0011-27-SP	10 x 27 x 11mm	0°
TSTS-MD0013-27-SP	10 x 27 x 13mm	0°
TSTS-MD0015-27-SP	10 x 27 x 15mm	0°
TSTS-MD0607-27-SP	10 x 27 x 7mm	6°
TSTS-MD0608-27-SP	10 x 27 x 8mm	6°
TSTS-MD0609-27-SP	10 x 27 x 9mm	6°
TSTS-MD0610-27-SP	10 x 27 x 10mm	6°
TSTS-MD0611-27-SP	10 x 27 x 11mm	6°
TSTS-MD0612-27-SP	10 x 27 x 12mm	6°
TSTS-MD0613-27-SP	10 x 27 x 13mm	6°
TSTS-MD0614-27-SP	10 x 27 x 14mm	6°
TSTS-MD0615-27-SP	10 x 27 x 15mm	6°
TSTS-MD0007-32-SP	10 x 32 x 7mm	0°
TSTS-MD0009-32-SP	10 x 32 x 9mm	0°
TSTS-MD0011-32-SP	10 x 32 x 11mm	0°
TSTS-MD0013-32-SP	10 x 32 x 13mm	0°
TSTS-MD0015-32-SP	10 x 32 x 15mm	0°
TSTS-MD0607-32-SP	10 x 32 x 7mm	6°
TSTS-MD0608-32-SP	10 x 32 x 8mm	6°
TSTS-MD0609-32-SP	10 x 32 x 9mm	6°
TSTS-MD0610-32-SP	10 x 32 x 10mm	6°
TSTS-MD0611-32-SP	10 x 32 x 11mm	6°
TSTS-MD0612-32-SP	10 x 32 x 12mm	6°
TSTS-MD0613-32-SP	10 x 32 x 13mm	6°
TSTS-MD0614-32-SP	10 x 32 x 14mm	6°
TSTS-MD0615-32-SP	10 x 32 x 15mm	6°