



CERVICAL SPINE TRUSS SYSTEM™ - STAND ALONE TITANIUM INTERBODY FUSION DEVICE

CSTS-SA has an advanced structural design that incorporates 4WEB Medical's proprietary Truss Implant Technology comprised of 3D printed titanium. Featuring integrated fixation and a quick anti-backout locking mechanism, the CSTS-SA is a spine solution for surgeon users in the stand-alone construct.



- ▲ Roughened titanium surface for an improved osteogenic response
- ▲ Web structure at the bone interface resists subsidence by distributing load across the entire vertebral endplate
- ▲ Open implant architecture allows for bone growth throughout the entire implant construct




Zero-profile stand-alone construct

Intuitive instrumentation available for guided and freehand surgical techniques

Quick anti-backout locking mechanism

Sterile packed for hospital efficiency and patient safety

5° of screw angulation




CSTS-SA IMPLANT SPECIFICATIONS

PRODUCT CATALOG

FOOTPRINT	HEIGHT	LORDOSIS
12 x 15mm	5-12mm	0°,7°
14 x 17mm	5-12mm	0°,7°

CATALOG #	FOOTPRINT	LORDOSIS
CSTS-SA-SM0005	12 x 15 x 5 mm	0°
CSTS-SA-SM0006	12 x 15 x 6 mm	0°
CSTS-SA-SM0007	12 x 15 x 7 mm	0°
CSTS-SA-SM0008	12 x 15 x 8 mm	0°
CSTS-SA-SM0009	12 x 15 x 9 mm	0°
CSTS-SA-SM0010	12 x 15 x 10 mm	0°
CSTS-SA-SM0011	12 x 15 x 11 mm	0°
CSTS-SA-SM0012	12 x 15 x 12 mm	0°
CSTS-SA-SM0705	12 x 15 x 5 mm	7°
CSTS-SA-SM0706	12 x 15 x 6 mm	7°
CSTS-SA-SM0707	12 x 15 x 7 mm	7°
CSTS-SA-SM0708	12 x 15 x 8 mm	7°
CSTS-SA-SM0709	12 x 15 x 9 mm	7°
CSTS-SA-SM0710	12 x 15 x 10 mm	7°
CSTS-SA-SM0711	12 x 15 x 11 mm	7°
CSTS-SA-SM0712	12 x 15 x 12 mm	7°

CSTS-SA-MD0005	14 x 17 x 5 mm	0°
CSTS-SA-MD0006	14 x 17 x 6 mm	0°
CSTS-SA-MD0007	14 x 17 x 7 mm	0°
CSTS-SA-MD0008	14 x 17 x 8 mm	0°
CSTS-SA-MD0009	14 x 17 x 9 mm	0°
CSTS-SA-MD0010	14 x 17 x 10 mm	0°
CSTS-SA-MD0011	14 x 17 x 11 mm	0°
CSTS-SA-MD0012	14 x 17 x 12 mm	0°
CSTS-SA-MD0705	14 x 17 x 5 mm	7°
CSTS-SA-MD0706	14 x 17 x 6 mm	7°
CSTS-SA-MD0707	14 x 17 x 7 mm	7°
CSTS-SA-MD0708	14 x 17 x 8 mm	7°
CSTS-SA-MD0709	14 x 17 x 9 mm	7°
CSTS-SA-MD0710	14 x 17 x 10 mm	7°
CSTS-SA-MD0711	14 x 17 x 11 mm	7°
CSTS-SA-MD0712	14 x 17 x 12 mm	7°



SELF-DRILLING SCREW SPECIFICATIONS

PRODUCT CATALOG

SCREW DIAMETER	LENGTH	
3.5mm	12,14,16mm	PRIMARY SCREW
3.8mm	12,14,16mm	RESCUE SCREW

CATALOG #	FOOTPRINT
CSCR-3512-SD	Ø3.5MM, L 12MM
CSCR-3514-SD	Ø3.5MM, L 14MM
CSCR-3516-SD	Ø3.5MM, L 16MM
CSCR-3812-SD	Ø3.8MM, L 12MM
CSCR-3814-SD	Ø3.8MM, L 14MM
CSCR-3816-SD	Ø3.8MM, L 16MM