




CERVICAL SPINE PLATING SOLUTION

4WEB's Cervical Spine Plating Solution (CSTS-PS) is compatible with the Cervical Spine Truss System™ which supports the cervical spine while biological fusion takes place. The CSTS-PS continues the legacy of delivering innovation, strength, and performance to every procedure.



- ▶ Wide variety of modular plating configurations to address multiple cervical spine pathologies and approaches.
- ▶ Plate design features a single-step locking mechanism to prevent screw backout.
- ▶ 1.8mm thick plates are available in 1-level, 2-level, 3-level, and 4-level configurations.




1-level



2-level



3-level




4-level




VISIT
PRODUCT
PAGE

CSTS-PS SPECIFICATIONS PLATES

CATALOG #	PLATE HEIGHT	HOLE-TO-HOLE HEIGHT
1-LEVEL PLATES		
CSTS-PLT116-NS	16mm	10mm
CSTS-PLT118-NS	18mm	12mm
CSTS-PLT120-NS	20mm	14mm
CSTS-PLT122-NS	22mm	16mm
CSTS-PLT124-NS	24mm	18mm
CSTS-PLT126-NS	26mm	20mm
CSTS-PLT128-NS	28mm	22mm
CSTS-PLT130-NS	30mm	24mm
CSTS-PLT132-NS	32mm	26mm
CSTS-PLT134-NS	34mm	28mm
CSTS-PLT136-NS	36mm	30mm
2-LEVEL PLATES		
CSTS-PLT228-NS	28mm	22mm
CSTS-PLT230-NS	30mm	24mm
CSTS-PLT232-NS	32mm	26mm
CSTS-PLT234-NS	34mm	28mm
CSTS-PLT236-NS	36mm	30mm
CSTS-PLT238-NS	38mm	32mm
CSTS-PLT240-NS	40mm	34mm
CSTS-PLT242-NS	42mm	36mm
CSTS-PLT244-NS	44mm	38mm
CSTS-PLT246-NS	46mm	40mm
CSTS-PLT248-NS	48mm	42mm
3-LEVEL PLATES		
CSTS-PLT342-NS	42mm	36mm
CSTS-PLT345-NS	45mm	39mm
CSTS-PLT348-NS	48mm	42mm
CSTS-PLT351-NS	51mm	45mm
CSTS-PLT354-NS	54mm	48mm
CSTS-PLT357-NS	57mm	51mm
CSTS-PLT360-NS	60mm	54mm
CSTS-PLT363-NS	63mm	57mm
CSTS-PLT366-NS	66mm	60mm
CSTS-PLT369-NS	69mm	63mm



CATALOG #	PLATE HEIGHT	HOLE-TO-HOLE HEIGHT
4-LEVEL PLATES		
CSTS-PLT463-NS	63mm	57mm
CSTS-PLT466-NS	66mm	60mm
CSTS-PLT469-NS	69mm	63mm
CSTS-PLT472-NS	72mm	66mm
CSTS-PLT475-NS	75mm	69mm
CSTS-PLT478-NS	78mm	72mm
CSTS-PLT481-NS	81mm	75mm



SELF-DRILLING SCREWS

SCREW DIAMETER	LENGTH	
3.5mm	12, 14, 16, 18mm	PRIMARY SCREW
3.8mm	12, 14, 16, 18mm	RESCUE SCREW
CATALOG #	DIAMETER, LENGTH*	
CSCR-3512-SD-SP	ø3.5mm, L 12mm	
CSCR-3514-SD-SP	ø3.5mm, L 14mm	
CSCR-3516-SD-SP	ø3.5mm, L 16mm	
CSCR-3518-SD-SP	ø3.5mm, L 18mm	
CSCR-3812-SD-SP	ø3.8mm, L 12mm	
CSCR-3814-SD-SP	ø3.8mm, L 14mm	
CSCR-3816-SD-SP	ø3.8mm, L 16mm	
CSCR-3818-SD-SP	ø3.8mm, L 18mm	