





## SI JOINT *TRUSS* SYSTEM™

4WEB's TRUSS Implant Technology™ redefined the standard of care as the first FDA-cleared 3D printed interbody fusion device. The SI Joint Truss System (SJTS) continues the legacy of delivering innovation, strength, and performance to every procedure.




- ▶ Dual pitch threads provide compression across the SI joint.
- ▶ 3D printed shell with titanium core to endure lateral and axial forces while fusion takes place.
- ▶ Cannulated and fenestrated with *TRUSS* structures to allow bone growth through the implant.
- ▶ Self-harvesting to maximize autograft utilization.


Instrumentation supports posterior and lateral approaches per surgeon preference



Threads extend to the tip allowing for ease of insertion




10.5, 12.0 and 13.5mm diameters in lengths from 35-90mm to accommodate patient anatomy



## PRODUCT CATALOG



VISIT  
PRODUCT  
PAGE



## SJTS IMPLANT SPECIFICATIONS

### THREADED SIJ FUSION DEVICE

CATALOG #	DIAMETER, LENGTH	CATALOG #	DIAMETER, LENGTH	CATALOG #	DIAMETER, LENGTH
SJTS-S-1035-SP	10.5mm x 35mm	SJTS-S-1235-SP	12.0mm x 35mm	SJTS-S-1335-SP	13.5mm x 35mm
SJTS-S-1040-SP	10.5mm x 40mm	SJTS-S-1240-SP	12.0mm x 40mm	SJTS-S-1340-SP	13.5mm x 40mm
SJTS-S-1045-SP	10.5mm x 45mm	SJTS-S-1245-SP	12.0mm x 45mm	SJTS-S-1345-SP	13.5mm x 45mm
SJTS-S-1050-SP	10.5mm x 50mm	SJTS-S-1250-SP	12.0mm x 50mm	SJTS-S-1350-SP	13.5mm x 50mm
SJTS-S-1055-SP	10.5mm x 55mm	SJTS-S-1255-SP	12.0mm x 55mm	SJTS-S-1355-SP	13.5mm x 55mm
SJTS-S-1060-SP	10.5mm x 60mm	SJTS-S-1260-SP	12.0mm x 60mm	SJTS-S-1360-SP	13.5mm x 60mm
SJTS-S-1065-SP	10.5mm x 65mm	SJTS-S-1265-SP	12.0mm x 65mm	SJTS-S-1365-SP	13.5mm x 65mm
* SJTS-S-1070-SP	10.5mm x 70mm	* SJTS-S-1270-SP	12.0mm x 70mm	* SJTS-S-1370-SP	13.5mm x 70mm
* SJTS-S-1075-SP	10.5mm x 75mm	* SJTS-S-1275-SP	12.0mm x 75mm	* SJTS-S-1375-SP	13.5mm x 75mm
* SJTS-S-1080-SP	10.5mm x 80mm	* SJTS-S-1280-SP	12.0mm x 80mm	* SJTS-S-1380-SP	13.5mm x 80mm
* SJTS-S-1085-SP	10.5mm x 85mm	* SJTS-S-1285-SP	12.0mm x 85mm	* SJTS-S-1385-SP	13.5mm x 85mm
* SJTS-S-1090-SP	10.5mm x 90mm	* SJTS-S-1290-SP	12.0mm x 90mm	* SJTS-S-1390-SP	13.5mm x 90mm

### THREADED SIJ FUSION DEVICE - LAG

CATALOG #	DIAMETER, LENGTH	CATALOG #	DIAMETER, LENGTH	CATALOG #	DIAMETER, LENGTH
* SJTS-L-1035-SP	10.5mm x 35mm	* SJTS-L-1235-SP	12.0mm x 35mm	* SJTS-L-1335-SP	13.5mm x 35mm
* SJTS-L-1040-SP	10.5mm x 40mm	* SJTS-L-1240-SP	12.0mm x 40mm	* SJTS-L-1340-SP	13.5mm x 40mm
* SJTS-L-1045-SP	10.5mm x 45mm	* SJTS-L-1245-SP	12.0mm x 45mm	* SJTS-L-1345-SP	13.5mm x 45mm
* SJTS-L-1050-SP	10.5mm x 50mm	* SJTS-L-1250-SP	12.0mm x 50mm	* SJTS-L-1350-SP	13.5mm x 50mm
* SJTS-L-1055-SP	10.5mm x 55mm	* SJTS-L-1255-SP	12.0mm x 55mm	* SJTS-L-1355-SP	13.5mm x 55mm
* SJTS-L-1060-SP	10.5mm x 60mm	* SJTS-L-1260-SP	12.0mm x 60mm	* SJTS-L-1360-SP	13.5mm x 60mm
* SJTS-L-1065-SP	10.5mm x 65mm	* SJTS-L-1265-SP	12.0mm x 65mm	* SJTS-L-1365-SP	13.5mm x 65mm
* SJTS-L-1070-SP	10.5mm x 70mm	* SJTS-L-1270-SP	12.0mm x 70mm	* SJTS-L-1370-SP	13.5mm x 70mm
* SJTS-L-1075-SP	10.5mm x 75mm	* SJTS-L-1275-SP	12.0mm x 75mm	* SJTS-L-1375-SP	13.5mm x 75mm
* SJTS-L-1080-SP	10.5mm x 80mm	* SJTS-L-1280-SP	12.0mm x 80mm	* SJTS-L-1380-SP	13.5mm x 80mm
* SJTS-L-1085-SP	10.5mm x 85mm	* SJTS-L-1285-SP	12.0mm x 85mm	* SJTS-L-1385-SP	13.5mm x 85mm
* SJTS-L-1090-SP	10.5mm x 90mm	* SJTS-L-1290-SP	12.0mm x 90mm	* SJTS-L-1390-SP	13.5mm x 90mm

### WASHERS\*\*

CATALOG #	IMPLANT SIZE
SJTS-W-1017	10.5mm
SJTS-W-1218	12.0mm
SJTS-W-1320	13.5mm

\* Available by request

\*\* Washer is sterile packaged with Threaded SIJ Fusion Device - Lag